




Oculus Orbit

SHAKING & MAGNETIC SEPARATION



Tailored for Lab Automation

Gentle orbital mixing: 1.5 mm radius orbital motion.
Adjustable speed of 0 – 1500 rpm for splatter-free mixing.

Plate locking: Automatic Plate Locking

Magnetic separation: Magnetic separation without moving the plate.
It automatically lifts the magnets under the plate.

Robust capacity: Supports loads of up to 1.000 g – from microplates to custom fixtures.

Customizable: Adapters for various panel sizes and magnetic configurations.

Let's build your solution together!

