



# Functionalized $\mu$ Beads

To capture your target

## Tailored for Lab Automation

**Gentle orbital mixing:** 1.5 mm radius orbital motion.

Adjustable speed of 0 – 1500 rpm for splatter-free mixing.

**Plate locking:** Automatic Plate Locking

**Magnetic separation:** Magnetic separation without moving the plate.  
It automatically lifts the magnets under the plate.

**Robust capacity:** Supports loads of up to 1.000 g – from microplates to custom fixtures.

**Customizable:** Adapters for various panel sizes and magnetic configurations.

**Let's build your solution together!**

